

Can You Put Any Amount of Embossing on a Tablet?

When it comes to tablet design, embossing is often viewed as an open canvas. Logos, brand names, dosage numbers, safety marks, or even decorative details can be applied to the punch face. This leads to a common misconception: that you can put any amount of embossing on a tablet face without running into problems. At first glance, the idea makes sense. If modern CAD software allows intricate designs, should they not translate directly into the tablet itself? However, the reality of manufacturing tablet tooling is different. Embossing is not limitless, and overstepping design boundaries can cause significant production and quality issues. Embossing sits at the intersection of engineering and material science, where choices have consequences that extend beyond aesthetics. In fact, even small changes to embossing depth or angle can have dramatic downstream effects on manufacturability and performance that are often overlooked at the concept stage.

Embossing has an important role in the veterinary, pharmaceutical, nutraceutical, and confectionery worlds. It is far more than cosmetic. By embossing a logo or symbol onto the surface of a tablet, manufacturers make products instantly recognisable to customers, pharmacists, and patients. This supports brand identity and consumer confidence, but it also provides critical functionality such as anti-counterfeiting, dosage identification, and even safety in preventing medication errors. In short, embossing is one of the most visible elements of tablet design, and it carries heavy responsibility. But while embossing is powerful, it is not without constraints. The geometry of the tablet, the behaviour of powders under compression, and the durability of the tooling all place limits on what is possible. Every embossing design therefore has to strike a careful balance between form and function, ensuring that aesthetic appeal never compromises mechanical integrity or patient safety.

One of the most fundamental limits is the available face area. The embossable portion of a tablet depends on its shape and profile. Flat-faced tablets and those with shallow cups offer generous surface area, sometimes up to 80 percent of the tablet face. However, more curved profiles, such as ball or pill shapes, reduce the usable area significantly, sometimes to only 40 percent. Attempting to emboss beyond these limits results in distorted characters, weakened tablet surfaces, or poor legibility. A flat logo may look perfect on a computer screen, but once applied to a small, curved tablet surface, the angles and depths of the strokes are altered, sometimes beyond recognition. This distortion becomes even more apparent when dealing with fine text or detailed logos, as perspective and curvature can subtly warp the perceived design. This distortion can become especially noticeable under certain lighting conditions, where even minor curvature can make embossed text appear uneven or incomplete. Ultimately, the flatter the tablet surface, the more freedom designers have to express visual identity without compromising manufacturability.

Excessive embossing also creates structural challenges. Every groove, cut, or recess in the tablet face represents a

stress point during compression. The more detail added, the weaker the tablet becomes. This can lead to manufacturing problems such as lamination, capping, or even tablet breakage during packaging and handling. It also affects the tooling itself. Deep, complex embossing can reduce the strength of the punch cup, making the punches more prone to wear or even catastrophic failure. In practice, embossing must be carefully balanced to preserve not just the integrity of the tablet, but also the longevity of the valuable tooling used to create it. Tooling engineers often describe this balance as a 'compromise of design intent,' where creative ambition must yield slightly to mechanical practicality. Even the most advanced tooling steels and coatings cannot compensate for fundamentally flawed embossing geometry, so good design remains the first line of defence against production failures. Manufacturers who adopt a proactive design-for-manufacture mindset often see fewer problems, faster setup times, and more consistent product quality.

Surface defects are another major concern when embossing is pushed too far. Problems such as picking and sticking are well-documented in tablet production. Picking occurs when powder particles cling to the fine details of an embossed design, pulling material out of the tablet face and leaving behind pits or voids. Sticking, on the other hand, is more general adhesion of product to the punch surface, often made worse by large, shallow embossing areas. Both issues compromise the appearance and quality of the tablet, and both become more likely as the density of embossing increases. In these cases, the problem is not the formulation alone, but the interaction between formulation, tooling, and the embossing design. Environmental factors such as humidity and temperature can exacerbate these problems. It highlights how even a seemingly minor design decision can have a ripple effect across the entire production line, from compression efficiency to final product presentation.

A good example of how embossing design directly influences tablet performance came from a customer who was struggling with persistent sticking. The issue appeared on both faces of the tablet, the upper face at the break-line and the lower face within the embossed detail. By working closely with the customer to review the design, it became clear that certain features were contributing to the problem. The tablet's soft centre was prone to adhesion, so a double radius profile was proposed, to remove the weak area and improve strength. The embossing itself contained trap zones that encouraged material build-up, so we simplified the design to eliminate those areas and improve release. To further reduce the risk of sticking, we introduced tapered peninsulas, which rounded off sharp, compound angles and allowed cleaner ejection. These design refinements improved the manufacture of the tablet, showing how practical design changes, developed in collaboration, resolve production issues and deliver functional and reliable tablets. The process also demonstrated the importance of cross-functional collaboration between tooling engineers, formulators, and production teams to achieve lasting results.

Even when tablets survive compression without major issues, film and sugar-coating processes can expose any weaknesses. If embossing is too shallow or wide, the coating



material may ‘infill,’ filling up the characters and erasing their definition. If the embossing is too steep, the coating may ‘bridge,’ over the grooves, leaving empty pockets beneath and producing weak or uneven surfaces. Both outcomes mean that the intended logo or text becomes blurred, illegible, or lost altogether. Worse still, coating stresses can erode the fine details, chip tablet edges, or even reduce tablet hardness. The idea that any embossing design can survive the rigors of coating simply does not hold up under real-world conditions. The coating process amplifies small flaws, so a design that appears perfect at compression may fail visually or mechanically once the coating has been applied and dried. Inconsistent coating coverage can also alter tablet appearance between batches, undermining brand uniformity and patient trust.

With embossed tool design, it is vital to recognise limitations and work within them. Successful tablet designs often start with simplifying logos or text to make them ‘tablet-friendly.’ This may mean opening up counters (the enclosed spaces inside letters like A, O, or P), avoiding sharp corners, or reducing the number of characters on a single face. Choosing the right tablet profile is also critical, since deeper cups reduce the available embossing area and distort strokes more severely. Specialised embossing fonts, developed specifically for pharmaceutical tooling, help minimise picking, sticking, and stress points while keeping characters clear and legible. Coating allowances must also be built into the design, with embossing depths adjusted to withstand the additional thickness of film or sugar layers. Even subtle adjustments, such as widening line spacing or rounding serif ends, can dramatically improve manufacturability without diminishing brand identity. These careful adjustments are often made during digital simulation phases, where virtual models can predict how the final tablet will behave under realistic conditions. This early attention to detail saves time, money, and frustration during later production stages.

When you consider all the individual factors, it becomes clear that embossing is already a highly technical design challenge. The complexity increases even further when multiple features are combined on a single tablet face. For example, pairing a company logo with dosage text and a functional break-line may seem straightforward on paper, but in practice each element competes for limited space


and interacts with the others in ways that affect strength, legibility, and manufacturability. A break-line might weaken the same area where embossing detail is placed, or a curved profile may distort one element more than another. The more design features added, the more careful the balance must be. Designers must also remember that patient usability is influenced by these factors; an unreadable dosage mark or an unclear logo could lead to confusion during dispensing or consumption. In regulated markets, such readability concerns are not just aesthetic but legal, since dosage misidentification can have safety implications. This makes the alignment of design, compliance, and manufacturing discipline more important than ever.

That is why it is always advised to work closely with tablet manufacturers at the very beginning of the design stage, reviewing every element in context and helping to identify the most effective combination of features. If this is not an option and manufacturing has already taken place using another design that has failed to meet expectations, we can still recommend targeted improvements that will achieve both brand impact and reliable manufacturability. Successful tablets are those where every element has been considered not just in isolation, but in agreement with tablet geometry, tooling strength, and downstream processes like coating. Collaboration early in the design process helps avoid costly tooling reworks, reduces time to market, and ensures that every batch maintains consistent visual and mechanical quality.

The myth that any amount of embossing can be placed on a tablet face ignores the realities of physics, tooling principles, and manufacturing. Overly ambitious designs risk producing weak, illegible, or inconsistent tablets, along with excessive tooling wear and coating problems. By contrast, carefully considered embossing, supported by collaboration with tooling experts, results in tablets that are strong, readable, and consistent across every production run. In essence, true design freedom comes not from ignoring limits but from understanding and respecting them.

The myth is busted. You cannot emboss without limits. Tablet design demands restraint and technical understanding. Far from being a blank canvas, the tablet face is a finely balanced surface where less often equals






PRODUCT NAME		MANUFACTURER		PRODUCT DESCRIPTION		COMPOSITION		TESTS		SPEC. & TESTS	
Varied Edge Embossing	ROUND	STANDARD	PEPPERMINT	1000 mg	1000 mg	CHOLESTEROL	0	SOA	2.21	HPGS 00.00N	HPGE 00.00N



more. The best embossing designs are those that achieve clarity, durability, and functionality without overloading the tablet with unnecessary detail. In tablet manufacturing, good design is not just about what you add, it is just as much about knowing where to stop. Designers who embrace these principles produce not only beautiful tablets but also efficient, reliable manufacturing outcomes that stand the test of time and market demands.



Andy Dumelow

Andy Dumelow joined I Holland as a PharmaCare® service engineer in 2000. He quickly moved into the role of trainer on the Approve Tooling Measuring System and was promoted to Product Manager, responsible for the running and development of the PharmaCare® range of products. He is now Head of Technical Sales for EMEA and NA.